


Titânio
Titanium
Ti 6AL 4v Eli
ISO 5832-3
ASTM F. 136



6212	CONTRA PARAFUSO Counter Screw
6214	CONECTOR TRANSVERSO Transverse Conector
6217	CONECTOR LATERAL Lateral Conector
6218	CONECTOR DE HASTE Rod Conector
6661	PLACA OCCIPITAL Occipital Plate
6668/01	GANCHO CURTO Short Hook
6668/02	GANCHO LONGO Long Hook

HASTE Ø 3,5 mm Rod	
REF.	Compr.mm Lenght
6663/01	80
/02	120
/03	240

PARAFUSO POLIAXIAL Ø3,5 mm Poliaxial Screw	
REF.	Compr. mm Lenght
6664/01	10
/02	12
/03	14
/04	16
/05	18
/06	20
/07	22
/08	24

PARAFUSO POLIAXIAL Ø 4 mm Poliaxial Screw	
REF.	Compr. mm Lenght
6665/01	10
/02	12
/03	14
/04	16
/05	18
/06	20
/07	22
/08	24
/09	26
/10	28
/11	30
/12	32
/13	34
/14	36
/15	38
/16	40
/17	42
/18	46
/19	50

PARAFUSO Ø 3,5 Screw		PARAFUSO Ø 4,0 Screw	
REF.	Compr.mm Lenght	REF.	Compr.mm Lenght
6662/01	6	6219/01	6
/02	8	/02	8
/03	10	/03	10
/04	12	/04	12
/05	14	/05	14
/06	16	/06	16

HASTE Ø 2,5 mm Rod	
REF.	Compr.mm Lenght
6213/01	30
/02	40

GANCHO CURVO Curved Hook	
REF.	Lado Side
6666/01	Direito Right
/02	Esquerdo Left