



**Inox**  
Stainless Steel  
ISO 5832-1  
ASTM F.138

**STEINMANN**

<i>LISO</i> <i>Plain</i>	<i>DIAMANTE</i> <i>Diamond</i>	COMPR. 230 mm <i>Length</i>		COMPR. 300 mm <i>Length</i>	
		REF.	Ø	REF.	Ø
		4680/01	3,5	4680/05	3,5
		/02	4,0	/06	4,0
		/03	4,5	/07	4,5
		/04	5,0	/08	5,0

<i>LISO</i> <i>Plain</i>	<i>TROCATER</i> <i>Trocar</i>	COMPR. 230 mm <i>Length</i>		COMPR. 300 mm <i>Length</i>	
		REF.	Ø	REF.	Ø
		4684/01	2,5	4684/05	2,5
		/02	3,5	/06	3,5
		/03	4,0	/07	4,0
		/04	5,0	/08	5,0

<i>ROSQUEADO</i> <i>Threaded</i>	<i>TROCATER</i> <i>Trocar</i>	COMPR. 230 mm <i>Length</i>		COMPR. 300 mm <i>Length</i>	
		REF.	Ø	REF.	Ø
		4688/01	2,5	4688/05	2,5
		/02	3,5	/06	3,5
		/03	4,0	/07	4,0
		/04	5,0	/08	5,0

**KIRSCHNER**

<i>LISO</i> <i>Plain</i>	<i>DIAMANTE</i> <i>Diamond</i>	COMPR. 230 mm <i>Length</i>		COMPR. 300 mm <i>Length</i>	
		REF.	Ø	REF.	Ø
		4728/01	1,0	4728/04	1,0
		/02	1,5	/05	1,5
		/03	2,0	/06	2,0
		/09	2,5	/10	2,5

<i>ROSQUEADO</i> <i>Threaded</i>	<i>DIAMANTE</i> <i>Diamond</i>	COMPR. 230 mm <i>Length</i>		COMPR. 300 mm <i>Length</i>	
		REF.	Ø	REF.	Ø
		4731/01	2,0	4731/02	2,0