


MINI FIN
Flexible Intramedullary Minimal Nail



Titânio
Titanium
Ti 6Al 7Nb
ASTM F. 1295

MINI FIN
Mini Nail

REF.	Ø mm	Tamanho	COR
Ref.		Size	Colour
5051/01	1,0	150mm	Azul - Blue
/02	1,2	150mm	Rosa- Pink
/03	1,5	150mm	Verde-Green



1734	CORTADOR MODELADOR DE FIOS <i>Cutter</i>
2171/01	CHAVE PARA INSERÇÃO DE FIO Ø1 A Ø1,5MM <i>T wrench with chuck</i>
2393	ALICATE CORTA FIO <i>Pliers</i>
2577	PERFURADOR <i>Perforator</i>