




ORTOSINTESE

GRAMPO E INSTRUMENTAL PARA EPÍFISE *Epiphyseal Staple And Instruments*



2177



2178



Inox

Stainless Steel
ISO 5832-1
ASTM F.138

GRAMPO TIPO BLOUNT Ø 2,4 mm
Staples Blount Type

REF.	LARGURA Width	ALTURA Heigth
4695	10	18
4696	10	22
4697	16	19
4698	16	22
4699	22	19
4700	22	22

2177 IMPACTOR INICIAL FINAL DE BLOUNT
Blount Initial Final Impactor

2178 EXTRATOR
Extractor

GRAMPO TIPO COVENTRY
Staples Coventry Type

REF.	LARGURA Width	ALTURA Heigth	mm
4701/01	24	33	5,0
/02	24	33	10,0
/03	24	33	15,0