




**ORTOSINTESE**

**FIXADOR EXTERNO PARA REDUÇÃO E ALONGAMENTO MODELO "ILIZAROV"**  
*External Fixator for Reduction and Alongation Model "Ilizarov"*




2261	SEMI ANEL DE FÊMUR 140° <i>Femur Semi Ring</i>
2262	SEMI ANEL DE FÊMUR 120° <i>Femur Semi Ring</i>
2263	SEMI ANEL DE FÊMUR 100° <i>Femur Semi Ring</i>
2264	SUPOORTE DE SEMI ANEL <i>Semi Ring Support</i>
2265	PORTA PINO DE SCHANZ PEQUENO <i>Small Schanz Pin Holder</i>
2266	PORTA PINO DE SCHANZ GRANDE <i>Large Schanz Pin Holder</i>

SEMI ANEL 5/8 <i>Semi Ring</i>	
REF.	Ø INTERNO Internal
2267/01	120
/02	130
/03	160

SEMI ANEL COM EXTREMIDADE ANGULADA <i>Semi Ring With Angled Extremity</i>	
REF.	Ø INTERNO Internal
2268/01	120
/02	130
/03	150

SEMI ANEL <i>Semi Ring</i>	
REF.	Ø INTERNO Internal
2260/01	100
/02	110
/03	120
/04	130
/05	140
/06	150
/07	160
/08	180
/09	200
/10	220
/11	240

**FIXADOR EXTERNO PARA REDUÇÃO E ALONGAMENTO MODELO "ILIZAROV"**  
*External Fixator for Reduction and Elongation Model "Ilizarov"*



2273	PLACA CURVA <i>Curved Plate</i>
2276	SUPORTE DE LIGAÇÃO COM ROSCA MACHO <i>Connection Support With Male Thread</i>
2277	SUPORTE DE LIGAÇÃO COM ROSCA FÊMEA <i>Connection Support With Female Thread</i>
2279	PARAFUSO FIXA FIO COM FURO PASSANTE <i>Wire Fix Screw With Passer-By Hole</i>
2280	PARAFUSO TENSOR DE FIOS COM CANAL <i>Wire Tightener Screw With Canal</i>
4753	PORCA SEXTAVADA M6 <i>Hexagonal Nut</i>

PLACA DE CONEXÃO ROSQUEADA <i>Threaded Connection Plate</i>		
REF.	FUROS Holes	COMPR. Length
2269/01	4	116
/02	7	176
/03	9	216
/04	11	256

PLACA DE CONEXÃO <i>Connection Plate</i>		
REF.	FUROS Holes	COMPR. Length
2270/01	8	155
/02	12	235
/03	17	335

BANDEIRA FÊMEA <i>Support With Female Thread</i>		
REF.	FUROS Holes	COMPR. Length
2274/01	2	40
/02	3	50
/03	4	60

BANDEIRA MACHO <i>Support With Male Thread</i>		
REF.	FUROS Holes	COMPR. Length
2275/01	2	27
/02	3	37
/03	4	47

PLACA DE CONEXÃO <i>Connection Plate</i>		
REF.	FUROS Holes	COMPR. Length
2271/01	3	45
/02	4	55
/03	5	65
/04	6	75
/05	7	85
/06	8	95
/07	9	105
/08	10	115

PLACA BI PLANA <i>Bi Plane Plate</i>		
REF.	FUROS Holes	COMPR. Length
2272/01	2	45
/02	3	65
/03	4	85

PARAFUSO COM CABEÇA SEXTAVADA  
*Hexagonal Screw With Head*

REF.	COMPR. Length
2278/01	15
/02	25
/03	30