



**ORTOSINTESE**

# SISTEMA ORTOLOCK FEMORAL CANULADO

## Ortolock Cannulated Femoral System

**Inox**  
Stainless Steel  
ISO 5832-1  
ASTM F 138

4841  
4842  
4843  
4844



3323



**PINO DE TRAVA**  
*Locking Bolt*

REF.	COMPR. Lengths
3178/05	34
/06	36
/07	38
/08	40
/09	42
/10	44
/11	46
/12	48
/13	50
/14	52
/15	54
/16	56
/17	58
/18	60
/19	64
/20	68
/21	72
/22	76
/23	80
/24	85
/25	90
/26	95
/27	100

3178



**ORTOLOCK FEMORAL CANULADO**  
*Cannulated Ortolock Femoral Nail*

	4841 Ø 10	4842 Ø 11	4843 Ø 12	4844 Ø 13	COMPR. Length
/01	●	●	●	●	300
/02	●	●	●	●	320
/03	●	●	●	●	340
/04	●	●	●	●	360
/05	●	●	●	●	380
/06	●	●	●	●	400
/07	●	●	●	●	420
/08	●	●	●	●	440
/09	●	●	●	●	460
/10	●	●	●	●	480

1221/02	GUIA ALINHADOR ROSQUEADO FEMORAL <i>Threaded Alignment Femoral Guide</i>
1222	GUIA APALPADOR <i>Alignment Guide</i>
2638/01	CABEÇA CONECTORA <i>Ortolock Guide</i>
3323	PARAFUSO TAMPÃO <i>Inox Plug Screw</i>

