



ORTOSINTESE


PLACA VOLFIX E PARAFUSO Volfix Plate and Screw

PLACA VOLFIX T OBLÍQUA ESQUERDA "T" Left Volfix Plate

REF.INOX	REF.Titânio	Furos
<i>Stainless Steel</i>	<i>Titanium</i>	<i>Holes</i>
6164/01	6137/01	3
/02	/02	4
/03	/03	6
/04	/04	8
/05	/05	10
/06	/06	12

PLACA VOLFIX T OBLÍQUA DIREITA "T" Right Volfix Plate

REF.INOX	REF.Titânio	Furos
<i>Stainless Steel</i>	<i>Titanium</i>	<i>Holes</i>
6163/01	6136/01	3
/02	/02	4
/03	/03	6
/04	/04	8
/05	/05	10
/06	/06	12



Inox
Stainless Steel
ISO 5832-1
ASTM F.139

Titânio
Titanium
Ti 6Al 4V Eli
ISO 5832-3
ASTM F. 136

PARAFUSO CORTICAL AUTO ROSQUEANTE Ø 3,5 mm Self Tapping Cortical Screw



REF. INOX	REF.TITÂNIO	COMPR.
<i>Stainless Steel</i>	<i>Titanium</i>	<i>Lenght</i>
4210/01	1361	10
/02	1362	12
/03	1363	14
/04	1364	16
/05	1365	18
/06	1366	20

PARAFUSO VOLFIX Ø2,5 ROSCA TOTAL Volfix Screw- Total Thread



INOX
Stainless Steel

Titânio
Titanium

6165 6151

PARAFUSO VOLFIX Ø2,5 ROSCA PARCIAL Volfix Screw Partial Thread



INOX
Stainless Steel

Titânio
Titanium

6166 6152

/01	12 mm
/02	14 mm
/03	16 mm
/04	18 mm
/05	20 mm
/06	22 mm
/07	24 mm
/08	26 mm
/09	28 mm
/10	30 mm
/11	32 mm
/12	34 mm
/13	36 mm
/14	38 mm
/15	40 mm
/16	42 mm
/17	44 mm
/18	46 mm
/19	48 mm
/20	50 mm
/21	52 mm
/22	54 mm
/23	56 mm
/24	58 mm
/25	60 mm

PARAFUSO CORTICAL AUTO ROSQUEANTE DE BLOQUEIO Ø 3,5 mm Cortical Screw Self Tapping With Blocking



REF. INOX	REF.TITÂNIO	COMPR.
<i>Stainless Steel</i>	<i>Titanium</i>	<i>Lenght</i>
4991	1400	10
4992	1401	12
4993	1402	14
4994	1403	16
4995	1404	18
4996	1405	20

PINO Ø2 Peg



REF.INOX	REF.Titânio	Tamanho
<i>Stainless Steel</i>	<i>Titanium</i>	<i>Size</i>
6167/01	6153/01	16 mm
/02	/02	18 mm
/03	/03	20 mm
/04	/04	22 mm
/05	/05	24 mm
/06	/06	26 mm
/07	/07	28 mm
/08	/08	30 mm