

PLACAS DE PEQUENOS FRAGMENTOS BLOQUEADAS

Small Fragments Locking Plates



Inox

Stainless Steel

ISO 5832-1

ASTM F.138

Inox

Stainless Steel

ISO 5832-1

ASTM F.138

Titânio

Titanium

Ti 6Al 4V Eli

ISO 5832-3


ASTM F. 136

Titânio

Titanium

ISO 5832-2

ASTM F. 67



PLACA DCP PEQ. FRAGMENTOS

Small Fragments DCP Plate

Titânio	AÇO INOX	FUROS
Titanium	Stainless Steel	Holes
6123/01	6570/01	4
/02	/02	5
/03	/03	6
/04	/04	7
/05	/05	8
/06	/06	9
/07	/07	10
/08	/08	11
/09	/09	12
/10	/10	14

PLACA EM "T"

"T" Plate

Titânio	AÇO INOX	FUROS
Titanium	Stainless Steel	Holes
6124	6571	3 x 3
6125	6572	3 x 5
6126	6573	3 x 7

Titânio	AÇO INOX	FUROS
Titanium	Stainless Steel	Holes
6127	6574	4 x 4
6128	6575	4 x 6

PLACA EM "T" OBLÍQUA DIREITA

Oblique "T" Plate (Right)

Titânio	AÇO INOX	FUROS
Titanium	Stainless Steel	Holes
6129/01	6576/01	3 x 3
/02	/02	3 x 4
/03	/03	3 x 5
/04	/04	3 x 7

PLACA EM "T" OBLÍQUA ESQUERDA

Oblique "T" Plate (Left)

Titânio	AÇO INOX	FUROS
Titanium	Stainless Steel	Holes
6130/01	6577/01	3 x 3
/02	/02	3 x 4
/03	/03	3 x 5
/04	/04	3 x 7

PLACA DE RECONSTRUÇÃO

Reconstruction Plate

Titânio	AÇO INOX	FUROS
Titanium	Stainless Steel	Holes
6131/01	6578/01	4
/02	/02	5
/03	/03	6
/04	/04	7
/05	/05	8
/06	/06	10
/07	/07	12

PLACA 1/3 TUBULAR

Tubular 1/3 Plate

Titânio	AÇO INOX	FUROS
Titanium	Stainless Steel	Holes
6133/01	4207/01	2
/02	/02	3
/03	/03	4
/04	/04	5
/05	/05	6
/06	/06	7
/07	/07	8
/08	/08	9
/09	/09	10
/10	/10	11
/11	/11	12