



ORTOSINTESE

PLACAS PEQUENOS FRAGMENTOS Plates Small Fragments

4420



4430



4446



4442



4439



Inox
Stainless Steel
ISO 5832-1
ASTM F.138

Inox
Stainless Steel
ISO 5832-1
ASTM F.139

4443



4438



PLACA "T" "T" Plate

REF.	FUROS Holes	COMPR. Length.
4438	4 + 4	56
4439	4 + 6	78
4442	3 + 3	50
4443	3 + 5	67

PLACA DE MÚLTIPLOS FRAGMENTOS Multiple Fragments Plate

REF.	FUROS Holes	COMPR. Length.
4446	7	45

• PARA USO COM PARAFUSO CORTICAL E ESPONJOSO DE 3,5 mm E PARAFUSO ESPONJOSO DE 4,0 mm

• For use With 3,5 mm Cortical and Cancellous Screw and 4,0 mm Cancellous Screw

PLACA AUTO - COMPRESSÃO Self - Compression Plate

REF.	FUROS Holes	COMPR. Length.
4420/01	2	25
/02	3	37
/03	4	49
/04	5	55
/05	6	73
/06	7	85
/07	8	79
/08	10	109
/09	12	121

PLACA DE 1/3 TUBO 1/3 Tube Plate

REF.	FUROS Holes	COMPR. Length.
4430/01	2	28
/02	3	40
/03	4	52
/04	5	64
/05	6	76
/06	7	88
/07	8	100
/08	9	112
/09	10	124