



Inox
Stainless Steel
ISO 5832-1
ASTM F.138

- Para uso com Parafuso Cortical de 4,5 mm e Parafuso Esponjoso de 6,5 mm.
- For use with 4,5 mm Cortical Screw and 6,5 mm Cancellous Screw

PLACA ESTREITA
Straight Plate

REF.	FUROS Holes	COMPR. mm Length
4003/01	03	54
/02	04	70
/03	05	86
/04	06	102
/05	07	118
/06	08	133
/07	09	149
/08	10	165
/09	12	197
/10	14	229
/11	16	260

PLACA LARGA
Broad Plate

REF.	FUROS Holes	COMPR. mm Length
4106/01	06	102
/02	07	118
/03	08	133
/04	09	149
/05	10	165
/06	12	197
/07	14	229
/08	16	260

PLACA SEMI TUBULAR 1/2 TUBO
Semi Tubular 1/2 Tube Plate

REF.	FUROS Holes	COMPR. mm Length
4202/01	02	39
/02	03	55
/03	04	71
/04	05	87
/05	06	103
/06	07	119
/07	08	135
/08	09	151
/09	10	167
/10	11	183
/11	12	199