



**Inox**  
Stainless Steel  
ISO 5832-1  
ASTM F.138

- PARA USO COM PARAFUSO CORTICAL DE 4,5 mm e PARAFUSO ESPONJOSO DE 6,5 mm
- For use With 4,5 mm Cortical Screw and 6,5 mm Cancellous Screw

PLACA ANGULADA PERTROCANTÉRICA 130°  
Pertrocanteric Angled Plate

REF.	FUROS Holes	PLACA Plate	LÂMINA Blade
4546/01	4	60	70
/02	4	60	80
/03	4	60	90
/04	6	104	60
/05	6	104	70
/06	6	104	80
/07	6	104	90
/08	9	150	60
/09	9	150	70
/10	9	150	80
/11	9	150	90
/12	12	200	60
/13	12	200	70
/14	12	200	80
/15	12	200	90

PLACA ANGULADA CONDILAR 95°  
Condylar Angled Plate

REF.	FUROS Holes	PLACA Plate	LÂMINA Blade
4534/01	5	90	60
/02	5	90	70
/03	5	90	80
/04	7	120	60
/05	7	120	70
/06	7	120	80
/07	9	150	60
/08	9	150	70
/09	9	150	80
/10	12	200	60
/11	12	200	70
/12	12	200	80