

PLACA TUBO E PARAFUSO DESLIZANTE SDQ e SDC Tube Plate And Sliding Screw DHS and DCS



PLACA TUBO 135° SDQ DHS Tube Plate

REF.	FUROS Holes
4571/01	3
/02	4
/03	5
/04	6
/05	7
/06	8
/07	9
/08	10
/09	12



Inox
Stainless Steel
ISO 5832-1
ASTM F.138

PLACA TUBO 95° SDC DCS Tube Plate

REF.	FUROS Holes
4589/01	6
/02	8
/03	10
/04	12

• PARA USO COM PARAFUSO CORTICAL DE 4,5 mm
E PARAFUSO ESPONJOSO DE 6,5 mm.

• Use with 4,5mm Cortical Screw and 6,5 mm Cancellous Screw.

PARAFUSO DESLIZANTE SDQ / SDC Sliding Screw For DHS / DCS Tube Plate

REF.	COMPR. mm	ROSCA mm
4594/01	51	20
/02	57	20
/03	63	20
/04	70	20
/05	76	20
/06	82	20
/07	89	20
/08	95	20
/09	102	20
/10	108	20
/11	114	20
/12	121	20
/13	127	20
/14	133	20
/15	140	20
/16	146	20
/17	152	20



PARAFUSO DESLIZANTE SDQ / SDC Sliding Screw For DHS / DCS Tube Plate

REF.	COMPR. mm	ROSCA mm
4595/01	51	30
/02	57	30
/03	63	30
/04	70	30
/05	76	30
/06	82	30
/07	89	30
/08	95	30
/09	102	30
/10	108	30
/11	114	30
/12	121	30
/13	127	30
/14	133	30
/15	140	30
/16	146	30
/17	152	30



ORTOSINTESE

PLACA TUBO BLOQUEADA E PARAFUSO DESLIZANTE *Locking Tube Plate And Sliding Screw*

PLACA TUBO 95° BLOQUEADA

DCS Locked Tube Plate

REF. FUROS Holes

4603/01	3
/02	4
/03	5
/04	6
/05	7
/06	8
/07	9
/08	10
/09	11
/10	12
/11	13
/12	14
/13	15
/14	16
/15	17
/16	18
/17	20
/18	22
/19	24



PLACA TUBO 135° BLOQUEADA

DHS Locked Tube Plate

REF. FUROS Holes

4604/01	2
/02	3
/03	4
/04	5
/05	6
/06	7
/07	8
/08	9
/09	10
/10	11
/11	12
/12	13
/13	14
/14	15
/15	16
/16	18
/17	20
/18	22
/19	24

Inox

Stainless Steel

ISO 5832-1

ASTM F.138

• PARA USO COM PARAFUSO CORTICAL DE 4,5 mm E PARAFUSO DE BLOQUEIO Ø 5,0,MM .PARA A PLACA REF.4589, UTILIZE PARAFUSO ESPONJOSO DE 6,5 mm.

• For use with Ø 4,5 mm Cortical Screw and Ø 5,0 mm Locking Screw .Use the Ø 6,5 mm Cancellous Screw for the plate Ref.4589.

PARAFUSO DESLIZANTE SDQ / SDC *Sliding Screw For DHS / DCS Tube Plate*

REF. COMPR. mm ROSCA mm

REF.	Compr. mm	ROSCA mm
4594/01	51	20
/02	57	20
/03	63	20
/04	70	20
/05	76	20
/06	82	20
/07	89	20
/08	95	20
/09	102	20
/10	108	20
/11	114	20
/12	121	20
/13	127	20
/14	133	20
/15	140	20
/16	146	20
/17	152	20

4594 4595



PARAFUSO DESLIZANTE SDQ / SDC *Sliding Screw For DHS / DCS Tube Plate*

REF. COMPR. mm ROSCA mm

REF.	Compr. mm	ROSCA mm
4595/01	51	30
/02	57	30
/03	63	30
/04	70	30
/05	76	30
/06	82	30
/07	89	30
/08	95	30
/09	102	30
/10	108	30
/11	114	30
/12	121	30
/13	127	30
/14	133	30
/15	140	30
/16	146	30
/17	152	30