


Inox
Stainless Steel
ISO 5832-1
ASTM F.139

Titânio
Titanium
Ti 6Al 4V Eli
ISO 5832-3
ASTM F.136



PLACA PARA ÚMERO LATERAL
Lateral Humerus Plate

REF.	REF	FUROS
Direita/Right	Esquerda/Left	Holes
6172/01	6173/01	3
/02	/02	5
/03	/03	7

PLACA PARA ÚMERO POSTERO LATERAL
Posterolateral Humerus Plate

REF.	REF	FUROS
Direita/Right	Esquerda/Left	Holes
6174/01	6175/01	3
/02	/02	5
/03	/03	7

PLACA OLÉCRANO
Olecranon Plate

REF.	REF	FUROS
Direita/Right	Esquerda/Left	Holes
6154/01	6155/01	2
/02	/02	3
/03	/03	4
/04	/04	5
/05	/05	6
/06	/06	7

PLACA PARA CLAVÍCULA DISTAL
Distal Clavicle Plate

REF.	REF	FUROS
Direita/Right	Esquerda/Left	Holes
6160/01	6161/01	4
/02	/02	5
/03	/03	6
/04	/04	7
/05	/05	8

PLACA PARA ÚMERO MEDIAL
Medial Humerus Plate

REF.	REF	FUROS
Direita/Right	Esquerda/Left	Holes
6170/01	6171/01	3
/02	/02	5
/03	/03	7

PLACA PARA ÚMERO POSTERO LATERAL C/ GANCHO
Posterolateral Humerus Plate with Hook

REF.	REF	FUROS
Direita/Right	Esquerda/Left	Holes
6176/01	6177/01	3
/02	/02	5
/03	/03	7

PLACA PARA CLAVÍCULA
Clavicle Plate

REF.	REF	FUROS
Direita/Right	Esquerda/Left	Holes
6157/01	6158/01	4
/02	/02	6
/03	/03	8

PLACA PARA CLAVÍCULA HOOK
Clavicle Hook Plate

REF.	REF	PROFUNDIDADE GANCHO
Direita/Right	Esquerda/Left	Depth Hook
6168/01	6169/01	12
/02	/02	15
/03	/03	18

PLACA PARA ÚMERO ANTERO MEDIAL
Anteromedial Humerus Plate

REF.	REF	FUROS
Direita/Right	Esquerda/Left	Holes
6178/01	6179/01	3
/02	/02	5
/03	/03	7

PLACA DE PUNHO
Wrist Plate

REF.	CURVATURA
6159/01	longa/long
/02	curta/ Short
/03	reta/Straight

PLACA SPIN
Spin Plate

Titânio	FUROS	COMP./mm	AÇO INOX	FUROS	COMP./mm
Titanium	Holes	Length	Stainless Steel	Holes	Length
6132/01	3	94	6580/01	3	94
/02	5	124	/02	5	124
/03	6	139	/03	6	139
/04	8	169	/04	8	169
/05	9	184	/05	10	199
/06	10	199	/06	12	229
/07	12	229	/07	9	184