


PRÓTESE DE QUADRIL TIPO CHARNLEY  
Hip Prosthesis Charnley Type



1015  
1016  
1017  
1018  
1019

1083

1028

**Inox**  
Stainless Steel  
ISO 5832-1  
ASTM F.138

1015/02	PRÓTESE MODULAR STANDARD 40 Modular Charnley Prosthesis
1015/03	PRÓTESE MODULAR STANDARD 40 POLIDA Modular Charnley Prosthesis (Polished)
1016/02	PRÓTESE MODULAR ESTREITA 35 Modular Charnley Prosthesis (Narrow)
1016/03	PRÓTESE MODULAR ESTREITA 35 POLIDA Modular Charnley Prosthesis (Polished Narrow)
1017/01	PRÓTESE 3/4 COLO MODULAR Modular 3/4 Neck Prosthesis
1018/01	PRÓTESE CDH Modular Modular CDH Prosthesis
1019/01	PRÓTESE EXTRA PEQUENA MODULAR Modular Charnley Prosthesis (XS)
1028	PRÓTESE CHARNLEY MODULAR Modular Flange Charnley Prosthesis

**CABEÇA INTERCAMBIÁVEL**  
*Interchangeable Head*

REF.	EXT. Ø mm	CAB Ø mm	Head Ext.
1083/01	22	MÉDIO	Medium
1083/02	22	LONGO	Long