




ORTOSINTESE

PRÓTESE DE QUADRIL SEM CIMENTO TIPO SPOTORNO *Cementless Hip Prosthesis Type Spotorno*



CABEÇA INTERCAMBIÁVEL *Interchangeable Head*

Cone STD	Ø mm	COMPR. COLO	Neck Length
1091/01	22	MÉDIO	Medium
/02	22	LONGO	Long
1088/01	28	CURTO	Short
/02	28	MÉDIO	Medium
/03	28	LONGO	Long
/04	28	EXTRA LONGO	Extra Long

COMPONENTE FEMORAL SEM CIMENTO *Cementless Femoral Component*

REF. TAMANHO/Size

1151/01	7
/02	8
/03	9
/04	10
/05	11,25
/06	12,50
/07	13,75
/08	15
/09	16,25
/10	17,5
/11	20

TAMANHO Size	A	OFFSET	α°
7	142,8	35	145°
8	146,4	36,1	145°
9	150	37,2	145°
10	153,6	38,2	145°
11,25	158,1	39,4	145°
12,5	162,6	40,6	145°
13,75	167,1	41,8	145°
15	171,6	43	145°
16,25	176,1	44,2	145°