

PRÓTESE DE QUADRIL SEM CIMENTO UNION

Cementless Hip Prosthesis Union



Titânio
Titanium
Ti 6AL 4V ELI
ISO 5832-3
ASTM F.136

1001



1088
1090
1091




Cr Co Mo
ISO 5832-4
ASTM F 75 ou
ASTM F 1537

CABEÇA INTERCAMBIÁVEL

Interchangeable Head

mm	Cone STD	CONE 12/14	COMPR.COLO	Neck Length
22	-	1091/05	CURTO	Short
22	1091/01	/06	MÉDIO	Medium
22	/02	/07	LONGO	Long
22		/08	EXTRA LONGO	Extra Long
28	1088/01	1088/05	CURTO	Short
28	/02	/06	MÉDIO	Medium
28	/03	/07	LONGO	Long
28	/04	/08	EXTRA LONGO	Extra Long
32	1090/01	1090/05	CURTO	Short
32	/02	/06	MÉDIO	Medium
32	/03	/07	LONGO	Long
32		/08	EXTRA LONGO	Extra Long

TAMANHO (SIZE)	"A"	"B"	"α"
9	130	41,3	135°
10	140	42,3	135°
11	145	42,8	135°
12	150	43,8	135°
13	155	44,3	135°
14	160	44,8	135°
15	165	45,8	135°
16	170	46,3	135°
18	180	47,3	135°



HASTE FEMORAL UNION

Union Femoral Stem

REF. TAMANHO

Cone STD	TAMANHO
1001/01	9
/02	10
/03	11
/04	12
/05	13
/06	14
/07	15
/08	16
/09	18

Cone 12/14	TAMANHO
/10	9
/11	10
/12	11
/13	12
/14	13
/15	14
/16	15
/17	16
/18	18