


PRÓTESE DE QUADRIL TIPO MULLER
Hip Prosthesis MullerType



Inox
 Stainless Steel
 ISO 5832-1
 ASTM F.138

PRÓTESE FEMORAL MODULAR
Modular Femoral Prosthesis

| REF. | Tamanho/Size |
|---------|--------------|
| 1080/01 | 7,5 |
| /02 | 10,0 |
| /03 | 12,5 |
| /04 | 15,0 |
| /05 | 17,5 |
| /06 | 20,0 |

CABEÇA INTERCAMIÁVEL
Interchangeable Head

| REF. | Ø mm | COMPR. COLO mm | Neck Length |
|---------|------|----------------|-------------|
| 1083/01 | 22 | NORMAL | Medium |
| | /02 | LONGO | Long |
| | /03 | EXTRA LONGO | Extra Long |
| 1082/01 | 28 | CURTO | Short |
| | /02 | NORMAL | Medium |
| | /03 | LONGO | Long |
| | /04 | EXTRA LONGO | Extra Long |
| 1081/01 | 32 | CURTO | Short |
| | /02 | NORMAL | Medium |
| | /03 | LONGO | Long |