





Cr Co Mo
ISO 5832-4
ASTM F.75




Titânio
Titanium
Ti 6AL 4V ELI
ISO 5832-3
ASTM F. 136

CABEÇA INTERCAMBIÁVEL
Interchangeable Head

| Cone STD | Ø mm | COMPR. | Length |
|----------|------|-------------|------------|
| 1091/01 | 22 | NORMAL | Medium |
| /02 | 22 | LONGO | Long |
| 1088/01 | 28 | CURTO | Short |
| /02 | 28 | NORMAL | Medium |
| /03 | 28 | LONGO | Long |
| /04 | 28 | EXTRA LONGO | Extra Long |

COLO PROXIMAL Ø 17
Proximal Component

| REF. | L3 | OFFSET |
|---------|----|--------|
| 1265/01 | 50 | 41 |
| /02 | 60 | 41 |
| /03 | 70 | 41 |
| /04 | 80 | 41 |
| /05 | 50 | 46 |
| /06 | 60 | 46 |
| /07 | 70 | 46 |
| /08 | 80 | 46 |



| | |
|------|--|
| 1270 | PARAFUSO TAMPÃO Ø10mm Plug Screw |
| 1271 | PARAFUSO DE FIXAÇÃO Ø 6 X 20mm Screws to fix the distal rod |

HASTE DISTAL
Distal Component

| REF. | Ømm | L1 | L2 |
|---------|-----|-------|-----|
| 1268/01 | 14 | 113,5 | 130 |
| /03 | 16 | 113,5 | 130 |
| /05 | 18 | 113,5 | 130 |
| /07 | 20 | 113,5 | 130 |
| /08 | 14 | 120 | 160 |
| /10 | 16 | 120 | 160 |
| /12 | 18 | 120 | 160 |
| /14 | 20 | 120 | 160 |
| /16 | 22 | 120 | 160 |
| /18 | 24 | 120 | 160 |
| /21 | 16 | 120 | 210 |
| /23 | 18 | 120 | 210 |
| /25 | 20 | 120 | 210 |
| /27 | 22 | 120 | 210 |
| /29 | 24 | 120 | 210 |

COLO PROXIMAL Ø 21.5
Proximal Component

| REF. | L3 | OFFSET |
|---------|----|--------|
| 1266/01 | 50 | 41 |
| /02 | 60 | 41 |
| /03 | 70 | 41 |
| /04 | 80 | 41 |
| /05 | 50 | 46 |
| /06 | 60 | 46 |
| /07 | 70 | 46 |
| /08 | 80 | 46 |